


# BIOQUELL SEQURE YOUR ULTIMATE SOLUTION FOR THOROUGH,

COMPLIANT DECONTAMINATION





Maintain your clean space's integrity with the Bioquell SeQure, a compact, fixed bio-decontamination system that provides increased productivity, robust regulatory compliance, and enhances cost and resource efficiencies through an automated, validated decontamination cycle using Hydrogen Peroxide Vapor.

#### EFFICACY IN A COMPACT WALL-MOUNTED SOLUTION



Using the latest advancements in automated decontamination technology, the Bioquell SeQure delivers complete decontamination of exposed surfaces via a wall-mounted system that is ideal for pass-throughs, material airlocks, and more. Bioquell's unique Hydrogen Peroxide Vapor technology reduces inherent risks, resources, and labor commonly associated with "spray and wipe" methods to decrease costs and potential contamination events introduced to your clean space.

#### FAST VALIDATION AND INSTALLATION



To uphold rigid safety requirements and GMP compliance without impeding lab or production efficiencies, the Bioquell SeQure system features automatic cycle prediction technology for fast validation. Whether installed in a brand-new facility or retrofitted into an existing workspace, the Bioquell SeQure can be installed and validated in facilities in less than one week. With decreased downtime and lower resource opportunity cost, your facility experiences a lessened financial impact and lower risk from potential contaminants.

### SMALL FOOTPRINT, LARGE IMPACT



The slim design and digital, user-friendly interface make the Bioquell SeQure ideal for use in many rooms and chambers, including connecting pass-throughs and material airlocks. Its minimal footprint is enabled by wall-mounting the components of the system inside or outside of the workspace, maximizing your working area for increased productivity.

bioquell.com

bic



# 

## YOUR BIOQUELL SEQURE SYSTEM IS OPTIMIZED WITH:

- · Bioquell SeQure Aeration Module
- · Equipment commissioning and IQ/OQ
- · Gas cycle development (GCD)
- Integration testing for Building Management System (BMS) and more complex systems
- Performance Qualification (PQ)
- User training
- · Complete documentation packages

#### **APPLICATIONS:**

You can regularly re-zero transferred equipment and supplies, as well as every surface in these areas:

- · Material airlocks and pass-throughs
- · Cleanrooms of varying levels
- Biotech / cell biology research labs
- $\cdot$  Bacterial / viral vector labs
- · Biopharma production facilities
- BSL-3 / BSL-4 labs

# 

## COMPONENTS:

#### Hydrogen Peroxide Vaporizer

- 4 integrated fans distribute vapor dynamically
- · Automated cycle prediction technology
- Operates at low temperatures, safe for heat-sensitive materials

#### **Control Unit**

- · Easy-to-use, color touchscreen
- One-touch technology initiates decontamination cycle



- Minimizes peroxide wastage with no decanting
- $\cdot$  Houses 2 H<sub>2</sub>O<sub>2</sub> bottles (950 ml each)

### GET TECHNICAL WITH BIOQUELL SEQURE

Physical and safety data	Operating data	Power data
Vaporizer measurements: 800 x 580 x 160 mm (31.5 x 22.8 x 6.3 in)	Temperature limit at start of cycles: 15-30°C (59-86°F) <sup>1</sup>	Regional power data requirements: 230v: single phase 50/60 Hz 6A 120v: single phase 50/60 Hz 15A 100v: single phase 50/60 Hz 15A
Bottle Module measurements: 520 x 550 x 140 mm (20.5 x 21.7 x 5.5 in)	Relative humidity (RH) limits at start of cycle: 80% max <sup>2</sup>	Operational consumption: 1.8 kW/h (max)
Wired Control Unit measurements: 260 x 350 x 95 mm (10.2 x 13.8 x3.7 in)	Hydrogen peroxide liquid 35% w/w Bioquell supply only <sup>3</sup> Fits 950 ml bottle (x2)	Supply: Installation category II

- The Bioquell SeQure will allow a cycle to start up to a maximum of  $35^{\circ}$  ( $95^{\circ}$ F) starting temperature. However, when operating above  $30^{\circ}$ C, please contact Bioquell to ensure the correct cycle is set up.
- <sup>2</sup> If requiring to operate outside this range, refer to Bioquell or its agents.
- <sup>3</sup> For a full hydrogen peroxide specification, please contact Bioquell or its agents.

#### The Americas

Bioquell Inc E: solutions@bioquell.com W: www.bioquell.com T: +1 215 682 0225

#### Global

UK Headquarters Bioquell UK Ltd **T:** +44 (0)1264 835 835

France Bioquell SAS T: +33 (0)1 43 78 15 94 **Germany** Bioquell GmbH **T:** +49 (0) 221 168 996 74

**Ireland** Bioquell Ireland **T:** +353 (0)61 603 622 China

Bioquell Shenzhen Technology Company Ltd **T:** +86 755 8635 2622

#### Singapore Bioquell Asia Pacific Pte Ltd T: +65 6592 5145



Disclaimer. Please note that this document comprises marketing literature and is for summary information purposes only; customers or potential customers must not rely upon the contents of this document. Bioquell, Inc. or its affiliates, distributors, agents or licensees (together 'Bioquell') reserve the right to make changes to the contents of this document at any time and without prior notification. Use Bioquell SeQure safely. Always read the label and product information before use. Bioquell is a registered trademark of Bioquell UK Ltd. © Bioquell, Inc. (2018). All rights reserved. LS001-MKT-156 US Rev 9