

# A leading laboratory testing company partners with Ecolab to save up to £1M per year by switching to Ecolab's Bioquell Qube

**Application:** Sterility Testing

**Segment:** Biopharmaceuticals

**Challenge:** Extensive downtime during and after bio-decontamination

**Recommended product:**  
Bioquell Qube



## Situation

**A world leader in laboratory testing was seeking an alternative to their current isolator, used for sterility testing, due to the following challenges:**

### ▼ Bio-decontamination Cycle

**Times:** On a monthly basis their isolator needed to be fully bio-decontaminated using hydrogen peroxide vapour. This bio-decontamination cycle, plus external exhaust, required a full day. Not being able to utilise the isolator created operational challenges, ultimately reducing production throughput.

### ▼ Employee Health and Safety

**Concerns:** The company sought to create a better working environment for its team members. Current isolators were cumbersome, uncomfortable, difficult to clean and required extended periods of standing while in a half suit.

### ▼ Operational Challenges & Cost:

The isolators resided in a Grade A cleanroom which created time delays for access and increased operational costs significantly. Additionally, it was important to reduce cost associated with cleaning, environmental monitoring and data analysis.

### ▼ Ongoing Maintenance Cost:

Maintaining the equipment also became a hassle as the company's staff performed all of the planned annual and unexpected maintenance on the isolators (excluding HEPA certifications). This had time consuming and costly consequences. Any alternative selected would need to provide maintenance and service work.



**TOTAL VALUE  
DELIVERED**

**£986,239**  
Savings Per Year

See page 3 for more details

## BIOQUELL QUBE

The Bioquell Qube is a modular isolator constructed of hardwearing polypropylene and integrated with Bioquell Hydrogen Peroxide Vapour bio-decontamination technology. The system comes in 2, 4 or 6 glove configurations and offers up to 2 material transfer chambers or rapid transfer ports (RTPs).



# A leading laboratory testing company partners with Ecolab to save up to £1M per year by switching to Ecolab's Bioquell Qube



## Solution

After a period of independent research, key decision makers met with Ecolab to review its Bioquell Qube, an isolator with integrated hydrogen peroxide vapour bio-decontamination technology.

The company subsequently purchased a Bioquell Qube system in an M-32 variation. This variation (pictured at the top of page 3) allows for 3 chambers of 2 gloves, and a material transfer hatch on each side of the system. This meant that up to two batches at once could be completed for sterility testing (2 chambers per system for testing, with one chamber used for bio-decontaminating any incoming supplies needed).

The Bioquell Qube's modular design concept allowed the company to select a unit that was customised for their operation. The Bioquell Qube was selected for the following reasons:

- ▾ Available in multiple configurations to suit their needs with up to 6 gloves and 2 transfer chambers
- ▾ Utilising integrated Bioquell hydrogen peroxide vapour bio-decontamination technology
- ▾ Operates in negative or positive pressure
- ▾ Validated for rapid cycles with a 6-log sporicidal kill on exposed surfaces (as short as 30 minutes)
- ▾ Optimised for sterility testing, including optional integrated sterility testing pump
- ▾ Software fitted to enable compliance with 21 CFR part 11 and Eudralex Vol. 4 Annex 11 regulations
- ▾ Does not require construction or electrical work for installation

## Results

The Bioquell Qube was built, installed and validated within 16 weeks from date of order. The installation, validation and training was managed by Ecolab's team of Installation Engineers and Validation Specialists.

They have been able to realise the following value:

- ▾ **Shorter Bio-decontamination Cycle Times:** Leveraging the Bioquell Qube's fully self-contained and validated hydrogen peroxide vapour bio-decontamination process, the total cycle time was reduced significantly from their previous isolator bio-decontamination process. The bio-decontamination step was reduced to only 3 hours, resulting in an additional 2-3 batches per week to be processed. The production value of this time was estimated at over £672,460 per year.
- ▾ **Improved Employee Health and Safety:** The Bioquell Qube has the user's comfort incorporated as part of the design allowing operators to stand or sit while working within the isolator chambers. No longer are half suits required due to the Bioquell Qube being a closed containment system. For the company, this meant that operators no longer needed to stand or wear a half suit for an extended period. Additionally, the Bioquell Qube sleeves are made from a lightweight, comfortable EPDM material ideal for extended use.

## ▾ Reduction in Operational

**Challenges and Cost:** With the Bioquell Qube able to be hosted in a lower grade cleanroom by creating a Grade A environment within its chamber(s), operational costs were able to be reduced and access to the isolator was much faster. As a result, they were able to save time and money. See how with these additional operational cost savings:

- ▾ Fewer supplies are used as they no longer require environmental monitoring after cycles, creating a supply saving of £10,410 per year.
- ▾ The Bioquell Qube has allowed for the company to use the Electronic Laboratory Notebook (ELN) for bio-decontamination cycles, saving time for data review. A daily reduction of 75 minutes per day analysing data results in savings of £121,040 per year.
- ▾ Lower operating costs for the lab space due to the Bioquell Qube's footprint compared to soft-walled isolators (and associated bio-decontamination units) creates a saving of £103,290 per year.
- ▾ When manual cleaning of the Bioquell Qube was required, it took a fraction of the time of the previous isolator, saving 2 hours per month or £6,420 per year.
- ▾ **Ongoing Maintenance Cost:** The client would also no longer need to perform maintenance on their equipment as the Bioquell Qube comes with a comprehensive service contract, saving £39,910 per year. Additionally, costs associated with HVAC integration were eliminated as the Bioquell Qube has a fully integrated HEPA process and does not require connecting to the building's HVAC system.



### Impact beyond cost savings

There are items that a dollar amount simply was not able to capture, but the impact was substantial:

- Shorter turnaround time (TAT) for high-risk samples has been a highly impactful benefit which allows them to better satisfy client requirements.
- A decrease in training time has made the staff more flexible, a benefit that can allow the company to meet client demand and improve client satisfaction as a result. The training period was significantly reduced by 50%.
- The time saved from using the ELN has been utilised on other revenue generating activities such as data review for client samples.

### Summary

With integrated and validated bio-decontamination using Bioquell Hydrogen Peroxide Vapour Technology, each bio-decontamination cycle would only make the isolator inaccessible for a fraction of the time compared to their previous isolators connected to an alternative bio-decontamination system that was not from Ecolab.

A validation process after installation ensured a 6-log sporicidal kill on every exposed surface within the chambers, creating an ideal, repeatable and monitored environment for sterility testing.

### ANNUAL SAVINGS

#### WASTE

**£10,410**  
saved on supply costs

#### PRODUCTIVITY

75 mins saved per day  
reviewing data equal to  
**£121,040**

2 hours saved per month on  
cleaning time equal to  
**£6,420**

**50%**  
decrease in  
training time

#### PROFITABILITY

2-3 extra sample batches  
per week equal to  
**£672,460**

**£39,910**  
saved on maintenance  
and services costs

#### ENERGY

**£103,290**  
saved on lower  
operating costs

#### TOTAL VALUE DELIVERED

**£986,239**  
in total value delivered per year

# How Ecolab adds value



## Site Surveys

Ecolab offers site-surveys from a trained technical team to review new or existing manufacturing and cleaning process.

### VALUE+

- Technical recommendations are matched to your most relevant objectives and key performance indicators (safety, water savings, efficiency, etc.)
- Formal report summarises potential savings
- Delivers implementation strategy that minimises production disruptions



## Lab Studies

Ecolab offers a laboratory team that can analyse soils, identify residues and test samples.

### VALUE+

- Creates a quicker and more successful cleaning recommendation
- Ensures the chemistry solution is effective against your unique, targeted soils
- Helps determine optimal dilution and implementation parameters



## Automated Bio-decontamination to meet your needs

Ecolab has automated bio-decontamination solutions to create a validated 6-log sporicidal kill using Bioquell Hydrogen Peroxide Vapour.

### VALUE+

- Save time and reduce waste by decreasing the risk of shutdown, production stoppage or product failure due to an environmental contaminant
- Help avoid time intensive investigations and potential patient and reputational harm

USE BIOQUELL PRODUCTS SAFELY. ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

#### ECOLAB LTD

52 Royce Cl  
Andover  
SP10 3TS, UK  
www.bioquell.com

#### EUROPE HEADQUARTERS

Richtistr. 7  
8304 Wallisellen  
Switzerland  
www.ecolab.com/lifesciences



## Validation Assistance

Ecolab can assist with the complex validation process and change controls by providing guidance and documentation.

### VALUE+

- Guidance includes:
  - Cleaning process design and development
  - Cleaning process validation
  - Continuous cleaning process verification
- Helps ensure process stays in line with regulatory expectations and acceptance criteria



## Dedicated Regulatory Team

Ecolab leverages industry expertise to deliver insights and guidance on regulatory expectations, trends and available tools to help ensure compliance for product safety and quality.

### VALUE+

- Helps ensure compliance with cGMP to ensure cleaning and disinfection meet product safety and quality requirements
- Help implement validations that meet regulatory expectations



## Customer Training

Ecolab's team is comprised of experts to help you implement, manage, and maintain quality & compliance requirements and protocols.

### VALUE+

- Gain valuable best-practices and implementation strategies from trained field and technical experts
- Reduce troubleshooting time from on-site support and process training to resolve problems quickly