

Bioquell HPV-AQ

Section: 1. PRODUCT AND COMPANY IDENTIFICATION

Product name : Bioquell HPV-AQ

Other means of identification Not applicable.

Recommended use Surface Disinfectant

Restrictions on use Reserved for industrial and professional use.

Product dilution information Product is sold ready to use.

Company Ecolab Ltd.

52 Royce Close, West Portway

SP10 3TS Andover United Kingdom

+44 (0) 1264 835 835, +44 (0) 1264 835 836

Bioquell.consumables@ecolab.com

Emergency telephone

number Issuing date : +1 760 476 3959. Use access code: 333809

05.05.2022

Section: 2. HAZARDS IDENTIFICATION

GHS Classification

Acute toxicity (Oral) : Category 4 Acute toxicity (Inhalation) : Category 4 Skin irritation Category 2 Serious eye damage Category 1

single exposure

Specific target organ toxicity - : Category 3 (Respiratory system)

GHS Label element

Hazard pictograms





Signal Word : Danger

Hazard Statements : Harmful if swallowed.

> Causes skin irritation. Harmful if inhaled.

Causes serious eye damage. May cause respiratory irritation.

Precautionary Statements : Prevention:

> Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray. Do not eat, drink or smoke when using this product. Wear protective gloves/ eye

protection/ face protection.

Response:

Immediately call a POISON CENTER or doctor/ physician. IF

921779

Version: 1.0 1/8 Revision Date :05.05.2022

Bioquell HPV-AQ

SWALLOWED: Call a POISON CENTER or doctor/ physician if you feel unwell. Rinse mouth. IF ON SKIN: Wash with plenty of soap and water. IF INHALED: Remove person to fresh air and keep comfortable for breathing. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Disposal:

Dispose of contents/container in accordance with local/regional/national/international regulations.

Other hazards : None known.

Section: 3. COMPOSITION/INFORMATION ON INGREDIENTS

Pure substance/mixture : Mixture

Chemical Name CAS-No. Concentration: (%)

Hydrogen peroxide 7722-84-1 30 - 60

Section: 4. FIRST AID MEASURES

In case of eye contact : Rinse immediately with plenty of water, also under the eyelids, for at

least 15 minutes. Remove contact lenses, if present and easy to do.

Continue rinsing. Get medical attention immediately.

In case of skin contact : Wash off immediately with plenty of water for at least 15 minutes. Use

a mild soap if available. Get medical attention if irritation develops and

persists.

If swallowed : Rinse mouth. Get medical attention if symptoms occur.

If inhaled : Remove to fresh air. Treat symptomatically. Get medical attention.

Protection of first-aiders : If potential for exposure exists refer to Section 8 for specific personal

protective equipment.

Notes to physician : Treat symptomatically.

Most important symptoms and effects, both acute and

delayed

: See Section 11 for more detailed information on health effects and

symptoms.

Section: 5. FIREFIGHTING MEASURES

Suitable extinguishing media : Water

Unsuitable extinguishing

media

: Anything other than water Carbon dioxide (CO2)

Dry chemical

Foam

Specific hazards during

firefighting

: Not flammable or combustible.

Hazardous combustion

products

: Depending on combustion properties, decomposition products may

include following materials:

Oxygen

921779

Version :1.0 2 / 8 Revision Date :05.05.2022

Bioquell HPV-AQ

for firefighters

Special protective equipment: Use personal protective equipment.

Specific extinguishing methods

: Fire residues and contaminated fire extinguishing water must be disposed of in accordance with local regulations. In the event of fire and/or explosion do not breathe fumes. Collect contaminated fire extinguishing water separately. This must not be discharged into drains.

Section: 6. ACCIDENTAL RELEASE MEASURES

Personal precautions. protective equipment and emergency procedures

: Ensure adequate ventilation. Keep people away from and upwind of spill/leak. Avoid inhalation, ingestion and contact with skin and eyes. When workers are facing concentrations above the exposure limit they must use appropriate certified respirators. Ensure clean-up is conducted by trained personnel only. Refer to protective measures listed in sections 7 and 8. Eliminate any possible source of ignition.

Environmental precautions

: Do not allow contact with soil, surface or ground water.

Methods and materials for containment and cleaning up : Stop leak if safe to do so. Contain spillage, and then collect with noncombustible absorbent material, (e.g. sand, earth, diatomaceous earth, vermiculite) and place in container for disposal according to local / national regulations (see section 13). Flush away traces with water. For large spills, dike spilled material or otherwise contain material to ensure runoff does not reach a waterway.

Section: 7. HANDLING AND STORAGE

Advice on safe handling

: Do not ingest. Do not get in eyes, on skin, or on clothing. Do not breathe dust/fume/gas/mist/vapours/spray. Use only with adequate ventilation. Wash hands thoroughly after handling. In case of mechanical malfunction, or if in contact with unknown dilution of product, wear full Personal Protective Equipment (PPE).

Conditions for safe storage

: Do not store on wooden pallets. Keep out of reach of children. Store in suitable labeled containers.

Storage temperature : 5 °C to 25 °C

Section: 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Components with workplace control parameters

Components	CAS-No.	Form of exposure	Permissible concentration	Basis
Hydrogen peroxide	7722-84-1	TWA	1 ppm 1.4 mg/m3	BH OEL
Hydrogen peroxide	7722-84-1	TWA	1 ppm 1.4 mg/m3	ARE OEL
Hydrogen peroxide	7722-84-1	TWA	1 ppm	ACGIH
		TWA	1 ppm 1.4 mg/m3	NIOSH REL
		TWA	1 ppm 1.4 mg/m3	OSHA Z-1

Bioquell HPV-AQ

Engineering measures : Effective exhaust ventilation system. Maintain air concentrations

below occupational exposure standards.

Personal protective equipment

Eye protection : Wear eye protection/face protection.

Hand protection : Recommended preventive skin protection

Gloves Nitrile rubber butyl-rubber

Breakthrough time: 1 – 4 hours

Gloves should be discarded and replaced if there is any indication of

degradation or chemical breakthrough.

Skin protection : No special protective equipment required.

Respiratory protection : When workers are facing concentrations above the exposure limit they

must use appropriate certified respirators.

Hygiene measures : Handle in accordance with good industrial hygiene and safety

practice. Wash face, hands and any exposed skin thoroughly after handling. Provide suitable facilities for quick drenching or flushing of

the eyes and body in case of contact or splash hazard.

Section: 9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid

Colour : clear, colourless

Odour : odourless

pH : 1.5 - 3.5, (100 %)

Flash point : Not applicable.

Odour Threshold : no data available

Melting point/freezing point : no data available

Initial boiling point and

boiling range

: > 100 °C

Evaporation rate : no data available
Flammability (solid, gas) : Not applicable.
Upper explosion limit : no data available
Lower explosion limit : no data available
Vapour pressure : no data available
Relative vapour density : no data available

Relative density : 1.1 - 1.2

Water solubility : soluble

Solubility in other solvents : no data available

Partition coefficient: n-

octanol/water

: log Pow: -1.57Method: Calculated

Version :1.0 4 / 8 Revision Date :05.05.2022

Bioquell HPV-AQ

Auto-ignition temperature : no data available

Thermal decomposition : no data available

Viscosity, kinematic : 0.980 mm2/s (40 °C)

Explosive properties : no data available

Oxidizing properties : Yes

Molecular weight : no data available VOC : no data available

Section: 10. STABILITY AND REACTIVITY

Reactivity : No dangerous reaction known under conditions of normal use.

Chemical stability : Contamination may result in dangerous pressure increases - closed

containers may rupture.

Possibility of hazardous

reactions

: No dangerous reaction known under conditions of normal use.

Conditions to avoid : Heat.

Exposure to light. Freezing temperatures.

Incompatible materials : Bases

Reducing agents Organic materials Strong oxidizing agents

Strong acids

Combustible material

Metals

Hazardous decomposition

products

: Depending on combustion properties, decomposition products may

include following materials:

Oxygen

Section: 11. TOXICOLOGICAL INFORMATION

Information on likely routes of : Inhalation, Eye contact, Skin contact

exposure

Potential Health Effects

Eyes : Causes serious eye damage.

Skin : Causes skin irritation.

Ingestion : Harmful if swallowed.

Inhalation : May cause respiratory tract irritation. Harmful if inhaled. May cause

nose, throat, and lung irritation.

Chronic Exposure : Health injuries are not known or expected under normal use.

Experience with human exposure

921779

Version :1.0 5 / 8 Revision Date :05.05.2022

Bioquell HPV-AQ

Eye contact : Redness, Pain, Corrosion

Skin contact : Redness, Irritation

Ingestion : No information available.

Inhalation : Respiratory irritation, Cough

Toxicity

Product

Acute oral toxicity : Acute toxicity estimate : 1,389 mg/kg
Acute inhalation toxicity : 4 h Acute toxicity estimate : 31.43 mg/l

Test atmosphere: vapour

Acute dermal toxicity : no data available
Skin corrosion/irritation : no data available
Serious eye damage/eye : no data available

irritation

Respiratory or skin

sensitization

: no data available

Carcinogenicity : no data available
Reproductive effects : no data available
Germ cell mutagenicity : no data available
Teratogenicity : no data available
STOT - single exposure : no data available
STOT - repeated exposure : no data available
Aspiration toxicity : no data available

Section: 12. ECOLOGICAL INFORMATION

Toxicity

Environmental Effects : This product has no known ecotoxicological effects.

Product

Toxicity to fish : no data available

Toxicity to daphnia and other : no data available

aquatic invertebrates

Toxicity to algae : no data available

Components

Toxicity to fish : Hydrogen peroxide

96 h LC50 Pimephales promelas (fathead minnow): 16.4 mg/l

Components

Toxicity to daphnia and other :

: Hydrogen peroxide

aquatic invertebrates

48 h LC50 Daphnia magna (Water flea): 2.4 mg/l

Components

Toxicity to algae : Hydrogen peroxide

921779

Version :1.0 6 / 8 Revision Date :05.05.2022

Bioquell HPV-AQ

72 h EC50 Skeletonema costatum (marine diatom): 1.38 mg/l

Persistence and degradability

Not applicable - inorganic

Bioaccumulative potential

no data available

Mobility in soil

no data available

Other adverse effects

no data available

Section: 13. DISPOSAL CONSIDERATIONS

Disposal methods : Do not contaminate ponds, waterways or ditches with chemical or

used container. Where possible recycling is preferred to disposal or

incineration. If recycling is not practicable, dispose of

contents/container in accordance with local regulations Dispose of

wastes in an approved waste disposal facility.

Disposal considerations : Dispose of as unused product. Empty containers should be taken to

an approved waste handling site for recycling or disposal. Do not reuse empty containers. Dispose of in accordance with local, state, and

federal regulations.

Section: 14. TRANSPORT INFORMATION

The shipper/consignor/sender is responsible to ensure that the packaging, labeling, and markings are in compliance with the selected mode of transport.

Land transport (ADR/ADN/RID)

UN number : 2014

UN proper shipping name : HYDROGEN PEROXIDE, AQUEOUS SOLUTION

Transport hazard class(es) : 5.1 (8)
Packing group : II
Environmental hazards : No
Special precautions for user : None

Air transport (IATA)

Not permitted for transport

Sea transport (IMDG/IMO)

UN number : 2014

Proper shipping name : HYDROGEN PEROXIDE, AQUEOUS SOLUTION

Class : 5.1 (8)
Packing group : II
Marine pollutant : No

Self-Accelerating : 60 °C

decomposition temperature

(SADT)

921779


Bioquell HPV-AQ

Section: 15. REGULATORY INFORMATION

For Regulatory information please contact your Ecolab sales representative.

Section: 16. OTHER INFORMATION

NFPA:



HMIS III:



0 = not significant, 1 = Slight, 2 = Moderate, 3 = High

4 = Extreme, * = Chronic

Issuing date : 05.05.2022

Version : 1.0

Prepared by : Regulatory Affairs

REVISED INFORMATION: Significant changes to regulatory or health information for this revision is indicated by a bar in the left-hand margin of the SDS.

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.